D

Gebrauchsinformation

LAVANOX-Serag® Wundspray

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SERAG-WIESSNER GmbH & Co. KG
Zum Kugelfang 8-12, D-95119 Naila

Hersteller:

Prontomed GmbH
Am Bahndamm 70
D-32120 Hiddénhausen
Made in Germany

GB

Instructions for use

LAVANOX-Serag® wound spray

Intended use:

LAVANOX-Serag® wound is used for cleansing, moisturising, mechanical decontamination and infection prophylaxis in acute and chronic wounds. During wound cleansing, LAVANOX-Serag® Wound Spray helps to remove biofilm and fibrinoid debris. In addition, LAVANOX-Serag® wound Spray supports a low-pain removal of wound dressings. Wound odors are quickly and reliably eliminated using LAVANOX-Serag® wound spray.

Product properties:

The application of LAVANOX-Serag® wound spray is non-irritating, and virtually painless. Residues of LAVANOX-Serag® wound Spray may remain in the wound. There is no restriction

on the frequency and duration of use. As a multi-dose container, LAVANOX-Serag® wound spray is preserved with NaOCl/HOCl. This allows the use of the opened individual product for 12 weeks.

Application information:

A wound should first be cleaned of dirt, blood and coatings. For this purpose, the affected areas are wetted with LAVANOX-Serag® wound Spray with several bursts of spray from a short distance and, if necessary, mechanically cleaned. In case of heavy contamination, the application can be repeated several times. After the application of LAVANOX-Serag® wound Spray, it is not necessary to rinse with another solution (i.e. residues of LAVANOX-Serag® wound Spray may remain in the wound). Encrusted and difficult to remove wound dressings or bandages are sprayed generously with LAVANOX-Serag® wound Spray until they are soaked. After a short waiting time, they can be carefully removed with little pain. If necessary, a thin layer of wound gel can then be applied. The wound gel and the wound spray gel from the LAVANOX-Serag® product line are suitable for this. LAVANOX-Serag® wound spray can be used without restriction regarding age and gender of the patient. There are no known undesirable side effects. A slight chlorine odor is a sign of the activity of the preservative system. LAVANOX-Serag® wound spray can be used by non-professionals if they have been instructed by healthcare professionals.

Application restrictions:

LAVANOX-Serag® wound spray should not be used if there are known allergies to any of the ingredients. For external use on skin wounds. Do not use in combination with other wound cleaning sprays, solutions, gels and pre-saturated dressings. Do not use simultaneously with products containing chlorhexidi-

ne. Do not use LAVANOX-Serag® wound spray on bleeding or exuding wounds. For use during pregnancy, while breast-feeding and in children under two years of age, consult your doctor or pharmacist beforehand. In children up to 12 years of age, the application should be supervised by an adult. If there is redness, swelling or pain in the wound area or if the condition of the wound does not improve within a few days, a doctor should be consulted.

Safety instructions:

Do not swallow. Avoid contact with eyes. Keep safe from children. Use only bottles that are in perfect condition and undamaged. Use the product only until the stated expiry date and not longer than 12 weeks after opening. Please report a potentially serious incident involving the product immediately to Prontomed (vigilanz@prontomed.de) and the competent authority in your country.

Ingredients:

Aqua, <0,08% sodium hypochlorite/hypochlorous acid

Issue: V01/2022.08

Distributed by:

SERAG-WIESSNER GmbH & Co. KG
Zum Kugelfang 8-12, D-95119 Naila

Manufacturer:

Prontomed GmbH
Am Bahndamm 70
D-32120 Hiddenhausen
Made in Germany

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