CE 1014

CE110203/N01925v3

de

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SERAPRO® RS, sterile
Surgical instrument made of
surgical steel

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both the appropriate techniques and the anatomy prior to use.
- The product is supplied sterile (sterilization method: ethylene oxide).
- The product is intended for single use only.
- Do not resterilize or re-use the product as critical changes in performance and safety are possible.
- Do not use after expiration date.
- Do not use if the product shows signs of damage (i.e. cuts, kinks, crushed area) or if outer package or sterile package are damaged or unintentionally opened before use.
- Do not bend the product into shape.
- Use the products only for the intended purpose.
- Please contact the manufacturer for any information.



STERILE EO

The attending physician is responsible for selecting the instruments for specific applications or rather surgical use.

INSTRUCTIONS FOR USE

It is important to ensure that parts of the textile implant are not completely pulled inside lumen of the tunneler SERAPRO® VT1. Non-observance may cause damage to the textile implant. Caution when handling the instrument. Damage to gloves/injuries are possible.

DESCRIPTION AND INTENDED PURPOSE

The product is a surgical instrument made of surgical stainless steel. It is used in combination with another surgical instrument (tunneler, SERAPRO® VT1; IS90) for the implantation of SERAG-WIESSNER textile implants in urogyneology.

C

possible. Surgical intervention may be necessary. Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS


The storage conditions apply as indicated on the label.

DISPOSAL

Disposal of unused products or waste material should be handled with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

D

SERAPRO® RS is guided through SERAPRO® VT1 (A). The arm of a textile implant is passed through the loop of SERAPRO® RS and fixed by moving the SERAPRO® RS slightly back (B and C). It is important that the arm of the implant is not completely pulled inside the lumen of the tunneler SERAPRO® VT1 (D).

CONTRAINDICATIONS

No product-related contraindications are known.

RESTRICTIONS AND LIMITATIONS OF USE

The products are not intended specifically for use in direct contact with the central nervous system, the heart or the central circulatory system.

UNDESIRABLE SIDE EFFECTS

Injuries to tissues, vessels, nerves or other body structures that can lead to bleeding and pain are

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SERAG-WIESSNER GmbH & Co. KG