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**Reprocessable instruments made of surgical steel****SERASIS®  
SERAPRO®****WARNINGS AND PRECAUTIONS**

- Read the instructions for use in their entirety before using the device.
- The device should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The device is supplied non-sterile and must be cleaned, disinfected and sterilised according to these instructions for use before the first and each subsequent use.
- The device must be checked for suitability for use and safety before each use.
- Do not use damaged device.
- Do not deform the device.
- Use the device only for the intended application.
- Please contact the manufacturer for further information.

**DESCRIPTION**

SERASIS® and SERAPRO® instruments are reusable instruments made of surgical steel. The attending physician is responsible for the selection of instruments for specific applications.

**INTENDED PURPOSE**

The medical device is classified as a surgical instrument used invasively during a surgical procedure. The instrument is intended for the implantation of textile implants in urology / urogynecology, which are also provided by SERAG-WIESSNER. Please refer to the separate product information.

**APPLICATION NOTE**

When selecting and using the surgical instrument, the condition of the patient, the experience of the professional personnel, the surgical technique as well as the type of the textile implant must be taken into account.

A safe handling of the devices is a prerequisite to prevent possible complications in association with them. This is why SERAG-WIESSNER offers workshops. Please contact us.

**CONTRAINDICATIONS**

There are no known contraindications.

**RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS**

This device is not intended for use in direct contact with the heart, the central circulatory system or the central nervous system.

The instructions listed under the item "Reprocessing" have been selected and validated by the manufacturer as a suitable reprocessing procedure for these medical devices. The

reprocessor is responsible for ensuring that the reprocessing actually carried out with the equipment, materials and personnel used in the reprocessing facility achieves the intended result. This requires validation and routine monitoring of the process.

### **UNDESIRABLE SIDE EFFECTS**

Injuries to tissues, organs, vessels, nerves or other body structures can occur, which can lead to bleeding and pain. Surgical intervention may become necessary.

Insufficient reprocessing can lead to wound contamination and / or wound infection.

Users and / or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

### **DISPOSAL INFORMATION**

Disposal of unused or damaged devices or waste material should be handled with particular care due to the potential risk of contamination and / or injury. In all other cases observe the specific national regulations.

### **REPROCESSING**

#### **General notes**

When selecting the reprocessing steps to be carried out, please take note of the details and regulations in this instructions for use as well as the national legal obligations and standards.

The devices affected by these instructions are to be classified in the risk level "Critical B" according to KRINKO / BfArM recommendations. Therefore, mechanical cleaning should always be carried out.

To ensure patient safety, the devices must be handled with care and protected from damage. Do not use devices that have been damaged in any way, as their performance and safety can no longer be guaranteed. Do not make any modifications or repairs.

Effective cleaning / disinfection is the indispensable prerequisite for effective sterilisation of the devices. Improper or inadequate cleaning can lead particles being deposited on the devices and reaching the area of intervention. In addition, functionality and lifetime of the device may be impaired.

The device packaging is not suitable for sterilisation!

The reprocessing should take place immediately. The instruments may only be reprocessed by qualified personnel who are able to assess the occurring risks with the corresponding effects.

Only validated procedures (device- and device-specific) may be used for cleaning / disinfection and sterilisation.

#### **Limitations on reuse**

The devices can, with appropriate care and provided they are undamaged and fully functional, be reprocessed and reused several times.

Surgical steel instruments prove to be robust and resistant to repeated reprocessing. The useful life is limited by:

- Functional limitation
- Loss of legibility of marking and coding
- Material degradation (hairline cracks, pitting, rust)

Surface discolorations are not a quality defect and do not lead to a functional limitation. The compatibility of the chemicals used with the devices must be guaranteed. Limitations can be found under "Material resistance".

### Instructions for reprocessing

<p><b>Initial treatment at the place of use</b></p>	<p>Coarse impurities should be removed immediately with a disposable cloth / paper tissue after each use. National requirements for hygiene measures should be observed.</p>
<p><b>Storage and transport</b></p>	<p>Storage and transport should be carried out with suitable containers, e. g. washing trays, to avoid damage to the devices and potential contamination and / or injury.</p>
<p><b>Pretreatment and preparation before cleaning</b></p>	<p>Remove all macroscopically visible contamination with a suitable brush (no metal brush or steel wool) under running water. Pay particular attention to the lumen to ensure that it is freely passable. If disassembly is necessary, it should be carried out according to the instructions. These can be found at the end of the instructions for use.</p> <ul style="list-style-type: none"> <li>• Do not use abrasive cleaning agents.</li> <li>• Do not place in saline solution (NaCl solutions).</li> <li>• There must be no damage caused by excessive force or sharp-edged objects.</li> <li>• Avoid overfilling instrument trays and wash trays.</li> <li>• The use of fixing agents or hot water (&gt; 40 °C) is prohibited.</li> <li>• Cleaning and disinfection should be carried out immediately after pretreatment (or within one hour at the latest).</li> </ul>
<p><b>Cleaning, disinfection and drying</b></p>	<p><u>Automatic cleaning and disinfection</u> The user is solely responsible for the correct compliance with all validated parameters.</p> <p>When selecting the cleaning agent to be used, please note the following:</p> <ul style="list-style-type: none"> <li>• The agent must be suitable for cleaning the device.</li> <li>• The chemicals used must be compatible with the device (for restrictions see "Material resistance").</li> </ul> <p>When selecting the disinfectant, the following should be noted:</p> <ul style="list-style-type: none"> <li>• It is recommended to use disinfectants that meet the requirements of the valid DIN EN ISO 15883 and have a proven effectiveness.</li> <li>• Only fully desalinated / demineralized water may be used for rinsing.</li> <li>• The disinfectant must be regularly maintained and checked according to internal and manufacturer's specifications.</li> </ul>

	<p>Follow the instructions of the detergent / disinfectant manufacturer. If alkaline cleaning agents are used, sufficient neutralisation must be carried out in accordance with the manufacturer's instructions.</p> <p><u>Automatic cleaning and disinfection process:</u></p> <ol style="list-style-type: none"> <li>1. Dismantle the dismountable instrument into individual parts according to the assembly instructions (see "assembly / disassembly instructions").</li> <li>2. Place the devices in the disinfectant basket or in a sterilisation tray. When fastening the devices, make sure that they do not touch each other. This impairs the effectiveness of the cleaning / disinfection.</li> <li>3. Place the disinfectant basket or sterilisation sieve with the devices in the disinfectant (Attention: When stacking several disinfectant baskets or sterilisation sieves, the manufacturer's instructions must be observed).</li> <li>4. Start the program: <ul style="list-style-type: none"> <li>• 1x pre-rinse with cold water without additives (2 min)</li> <li>• Cleaning at 55 °C for at least 5 min with a suitable cleaning agent</li> <li>• Intermediate rinsing with cold water and neutralisation, if alkaline cleaner was used</li> <li>• Intermediate rinsing with cold water</li> <li>• Thermal disinfection at 90 – 93 °C for at least 5 min (A<sub>0</sub>-value &gt; 3000)</li> </ul> </li> <li>5. Remove the disinfectant basket or the sterilisation tray with the devices after the end of the programme. When removing the devices, make sure that any visible impurity is completely removed, especially in recesses, cavities, lumens and openings. If necessary, repeat the operation or use manual cleaning. Clean lumens and channels with soft round plastic brushes. The diameters of the lumen and brush must match.</li> <li>6. In case of residual moisture - place devices in drying cabinet at 50 °C to 60 °C for subsequent drying. The drying time depends on the load and the type of dishes.</li> <li>7. Pack the devices as soon as possible after removing them from the disinfectant or after maintenance, inspection and testing (see "Packaging" instructions).</li> </ol>
<p><b>Maintenance, inspection and testing</b></p>	<p>After each cleaning / disinfection, the devices must be macroscopically clean (including in all cavities), i.e. free of visible protein residues and other contaminants. In case of residual contamination, the complete cleaning and disinfection process must be repeated.</p> <p>Worn, damaged and porous devices must be sorted out.</p> <p>Corroded devices must also be removed, as they can cause corrosion on intact devices through extraneous rust transfer.</p>

	<p>The medical device must be assembled according to the assembly instructions and its functionality must be checked.</p> <p>Instrument oils should not be used if possible. If use is nevertheless desired, care should be taken to ensure that only instrument oils (Liquid paraffin) are used which - taking into account the maximum sterilisation temperature applied - are approved for steam sterilisation and have a tested biocompatibility.</p> <p>Pack the devices immediately (see "Packaging" instructions).</p>
<b>Packaging</b>	<p>Before sterilisation, the devices must be placed in suitable packaging. The packaging must be suitable for steam sterilisation (temperature resistance, sufficient steam permeability).</p> <p><u>Single packed:</u> A standardised packaging material can be used. The packaging must be large enough so that a seal is not under tension.</p> <p><u>Sets:</u> Sort devices into trays provided or place them on sterilisation trays. Cutting edges should be protected. A suitable procedure should be used for packing the trays.</p>
<b>Sterilisation</b>	<p>The specified times and temperatures are minimum requirements which must not be undercut. Should a downward deviation be necessary for process-technical reasons, this is the sole responsibility of the user and must be validated by the user.</p> <p>An exceeding of the specified times and temperatures is generally possible. However, longer sterilisation times and higher temperatures lead to increased stress on the material, which can result in premature aging of the devices.</p> <p>Only cleaned and, if necessary, disinfected devices may be sterilised. When using steam sterilisers, it must be ensured that sterilisation steam without contamination is used.</p> <p>Sterilisation is carried out with a steam steriliser, which has a fractionated vacuum process. Other sterilisation methods must be validated by the user.</p> <p>The exposure time is at least 3 min at a sterilisation temperature of 132 °C. The maximum sterilisation conditions are 18 min and 135 °C.</p>
<b>Storage</b>	<p>Sterilised devices should be stored in a dry and clean environment. The storage time depends on the type of packaging. Information on this can be found in the recommendation on the storage time for sterile medical devices of the Technical Committee Quality of the DGSV or evaluated accordingly.</p>
<b>Additional Information</b>	<p>If the specified chemicals and machines are not available, it is the responsibility of the user to validate his process accordingly.</p> <p>The user must ensure that the reprocessing process (including resources, materials and personnel) is suitable to achieve the required results. The state of the art and national laws require that</p>


	<p>validated processes are followed. The validated parameters must be maintained during each reprocessing cycle. In addition, the respective valid national legal regulations as well as the hygiene regulations of the respective health facility apply.</p> <p>If devices are used on patients suspected of having a prion disease, there is a potentially high risk of transmission after surgery. Whether the device must be destroyed or can be reprocessed is entirely at the discretion of the user.</p> <p>Defective devices that are returned to the manufacturer must first have gone through the entire recycling process.</p>
<b>Material resistance</b>	<p>When selecting cleaning agents and disinfectants, make sure that they <b>do not</b> contain the following ingredients:</p> <ul style="list-style-type: none"> <li>• organic, mineral and oxidising acids</li> <li>• oxidising agents (e.g. hydrogen peroxide).</li> </ul>
<b>Materials, equipment and parameters used during validation</b>	<p><b>Automatic cleaning and disinfection</b></p> <p>Proof of the basic suitability of the devices for effective automatic cleaning and disinfection has been provided by an independent test laboratory using the following materials, equipment and parameters:</p> <p><u>Washer-disinfector:</u> Miele Professional G 7836 CD, Miele &amp; Cie. KG, Gütersloh.</p> <p><u>Method:</u> Pre-cleaning with cold water 2 min; cleaning, 55 °C for 5 min, detergent concentration 0.5 %; rinsing with cold water 5 min; disinfection: 90 °C for 5 min.</p> <p><u>Cleaning agent:</u> <i>neodisher® MediClean forte</i>, Dr. Weigert GmbH &amp; Co KG, Hamburg, Ref 405033.</p> <p>The procedure described under the instruction "Cleaning and disinfection" (automatic) was followed.</p> <p><b>Packaging</b></p> <p>The following materials and equipment have been used in accordance with the procedure described under "Packaging": STERIKING® See-through reel, flat roll, type R43 and R44, Wipak.</p> <p><b>Sterilisation</b></p> <p>Proof of the basic suitability of the devices for effective steam sterilisation has been provided by an independent test laboratory using the following materials, equipment and parameters:</p> <p><u>Steam steriliser:</u> <i>Selectomat HP 666-1 HRED (serial no. SN 971118)</i>, MMM Group, Planegg.</p> <p><u>Method:</u> fractionated vacuum process.</p> <p>Exposure time: at least 3 min. Sterilisation temperature: 132 °C.</p> <p>Typical conditions in clinics and medical practices were taken into account and the procedure described under the instruction "Sterilisation" was followed.</p>

<b>Validation</b>	<p>The reprocessing process described in this document (pre-treatment, automatic and manual cleaning / disinfection, sterilisation) has been validated by an independent testing laboratory. Details of the validated process can be found under "Materials, equipment and parameters used in the validation".</p> <p>For the instructions mentioned here, the following standards were considered by the testing laboratory during validation:</p> <p>Cleaning: EN ISO 15883-1 Disinfection: EN ISO 15883-1 and -2 Packaging: EN ISO 11607-1 Sterilization: DIN EN ISO 14937, DIN EN 285 and DIN EN ISO 17665-1</p> <p>The validation was carried out using the following devices: SERAPRO® VT1 (REF: IS90), clip applicator (DV2010-01). With reference to the group or family validation mentioned in DIN EN ISO 17664, the present reprocessing process can be applied to all reusable surgical devices from SERAG-WIESSNER - under consideration of the worst case in connection with the proof of common characteristics and attributes.</p>
-------------------	--

## MONTAGE- / DEMONTAGEANLEITUNG / ASSEMBLY / DISASSEMBLY INSTRUCTIONS

### IS23 – Montage des Instruments / Assembly of the medical device

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## IS23 – Demontage des Instruments / Disassembly of medical device

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## IS74 – Montage des Instruments / Assembly of the medical device

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## IS74 – Demontage des Instruments / Disassembly of the medical device

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## IS200 – Montage des Instruments / Assembly of the medical device

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### IS200 – Demontage des Instrumentes / Disassembly of the medical device

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## IS210 – Montage des Instruments / Assembly of the medical device

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## IS210 – Demontage des Instruments / Disassembly of the medical device

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## IS211 – Montage des Instruments / Assembly of the medical device

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## IS211 – Demontage des Instruments / Disassembly of the medical device

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SERAG-WIESSNER GmbH & Co. KG