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DESCRIPTION

LAVANID® 1 and LAVANID® 2 are sterile, isotonic wound irrigation solutions based on Ringer's solution, containing 0.02% and 0.04% polyhexanide

respectively as preservative. Polyhexanide prevents the growth of micro-organisms in wound dressings and in the solution.

COMPOSITION

Ringer's solution (isotonic electrolyte solution made up of sodium chloride, potassium chloride, calcium chloride 2 H₂O, water for injection), polyethylene glycol, polyhexanide 0.02% (LAVANID® 1) or polyhexanide 0.04% (LAVANID® 2).

INDICATIONS

- Solution for local wound irrigation and cleansing
- To moisten cloths, tamponades and dressings

CONTRAINDICATIONS

LAVANID® must not be used:

- on any part of the CNS (including the meninges)
- in the area of cartilage and joints
- in the middle or inner ear
- intraperitoneally
- in cases of known hypersensitivity to any of the ingredients
- during pregnancy and breast-feeding

WARNINGS AND PRECAUTIONS FOR USE

Use clear solutions only!

Do not mix with additives.

Do not resterilize! This may critically alter its usability. The stoppered closure allows the aseptic removal of irrigation solution, e.g. using a spike top.

Alternatively, the screw cap can be removed so that the solution can be poured out.

INTERACTIONS WITH OTHER DRUGS

When used in combination with other wound treatment products, the instructions for use of these products should be checked for possible interactions with the ingredients of LAVANID®.

Due to incompatibility problems with the ingredient polyhexanide, LAVANID® must not be used in combination with anionic substances (e.g. anionic tensides or anionic gel formers).

Preoperative, surface skin disinfection with PVP-iodine preparations and subsequent intraoperative irrigation with LAVANID® is possible. Avoid direct contact between PVP-iodine preparations and the LAVANID® irrigation solution in the same wound area.

SIDE EFFECTS

The ingredient polyhexanide can cause allergic reactions. There are individual case reports of acute immune system reactions (anaphylactic reactions). Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

POSOLOGY, METHOD AND DURATION OF ADMINISTRATION

For wound irrigation during surgical procedures and diagnostic and therapeutic measures that require irrigation or moistening.

INFORMATION ABOUT SHELF-LIFE

The expiry date can be found on the label. Do not use after the expiry date!

Shelf life after opening:
2 months.

DISPOSAL INFORMATION

Dispose of unused products or waste material in accordance with the specific national regulations.

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SERAG-WIESSNER GmbH & Co. KG