# CE 1014

CE110105 / N01924v5

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### INSTRUCTIONS FOR USE

### WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified personnel, that is responsible to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).
- The product is intended for single use only.
- Do not resterilize or re-use the product as critical changes in performance and safety are possible.
- Do not use after expiration date.

- Do not use if the product shows signs of damage (i.e. cuts, kinks, crushed area) or if outer package or sterile package are damaged or unintentionally opened before use.

- Do not bend the product into shape.
- Use the products only for the intended purpose.

## DESCRIPTION AND INTENDED PURPOSE

The products are surgical eye needles made of stainless steel. The products are offered in different variants with regard to shape, length and tip. Eye needles serve as surgical suture material in combination with sterile thread material.

The attending physician is responsible for selecting the instruments for specific applications or rather surgical use.

## INSTRUCTIONS FOR USE

The thread material is inserted into the eye of the needle.

In eye needles the thread is passed with the long side backwards over the insertion funnel of the upper eye through the slit opening, whereby the needle should be clamped in a needle holder. A sufficient extension of the thread prevents unintentional removal. If a higher degree of security is required, the thread can be secured against slipping out by threading it into the lower eyelet.

The thread assignment table below shows the recommended needle-thread combinations. Non-observance may cause damage to the thread or needle. Needlestick injuries are possible.

To avoid needle damage during suturing, the needle should be grasped from the end of the needle (eye side) in an area between the end of the first third and the middle of the needle seen from the end of the needle (eye side).

The deformation of needles or holding them outside the recommended range may result in result in damage or needle breakage.

THREAD ASSIGNMENT		
Gauge size USP	Needle type	Gauge size EP
5/0 – 4/0	DSF-16	1 – 1.5
5/0 – 3/0	DS-14 – DS-33 DSF-21 – DSF-36 GR-22 – GR-36 HR-15 – HR-22 HR-15-F – HR-44-F HRF-20 – HRF-22 HS-15 – HS-20 HSF-12 – HSF-17	1 – 2
2/0 – 0	DS-36 – DS-40 GR-60 HR-26 – HR-38 HRF-26 – HRF-40 HRT-20 – HRT-26 HRX-20 HS-23 – HS-33 HSF-20 HSX-20	3 – 3.5
2/0 – 1	HRF-45 – HRF-70	3 – 4
0 – 2	DS-43 – DS-63 GS-42 – GS-45 HR-42 – HR-55 HRT-30 – HRT-38 HRX-25 – HRX-40 HS-36 – HS-55 HSX-25 – HSX-38	3.5 – 5
1 – 3+4	DS-68 – DS-95 HR-62 – HR-95 HRT-42 – HRT-85 HRX-46 – HRX-55 HS-60 – HS-75 HSX-40 – HSX-48	4 – 6
2 – 5	HRX-68 HSX-55 – HSX-70	5 – 7
3+4 – 8	HSX-75	6 – 10

## CONTRAINDICATIONS

No product-related contraindications are known.

## RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The products are not intended specifically for use in direct contact with the central nervous system, the heart or the central circulatory system.

## ADVERSE REACTIONS

Injuries to tissues, vessels, nerves or other body structures that can lead to bleeding and pain are possible. Surgical intervention may be necessary. Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

## STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

## DISPOSAL INFORMATION

Disposal of unused products or waste material should be handled with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

## MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG